

**SOTER™ CONDENSER (#8777)**



Note: Overall System Set-Up Area Requirement with optional SOTER™ Condenser is 36" (W) x 30" (D) 92cm (W) x 76cm (D)



**# 8777 SOTER™ Condenser**

SOTER™ Condenser quietly separates the waste 'Purge' alcohol and CO<sub>2</sub> discharge from your SAMDRI® process chamber. Comes complete with 15' static free flexi-exhaust hose w/ 1/4" ID tubing, an insulated exhaust connection from SAMDRI® to it, and clear drain flex tubing. Both 120V and 220V availability.

**Benefits**

1. Safely separates the waste Alcohol and CO<sub>2</sub> in a safe manner.
2. Eliminates static discharge that may occur during the PURGE mode of the process.
3. 100% noise cancelling and containment of any exhaust discharge during PURGE mode before safely routing to Carboy.



**SOTER™ Condenser Shown Connected to SAMDRI®-795**