

---

 **tousimis** Large Capacity Holders: 4", 3", 2" and 1"



... Excellent for Gel, Large Biological or Material Samples ...

# Large Capacity Holders



**4" Large Capacity Holder (Material: SS316) - Cat# 8771D**  
Interior 0.7" Ht x 4.0" I.D. (17.8 mm Ht x 101.6 mm I.D.)  
Exterior 1.05" Ht x 4.4" O.D. (26.7 mm Ht x 111.8 mm O.D.)  
Mesh Material: SS 316 / Mesh Size: 35x35 / Opening: 0.017" (0.43mm)



**3" Large Capacity Holder (Material: SS316) - Cat# 8771C**  
Interior 0.7" Ht x 3.0" I.D. (17.8 mm Ht x 76.2 mm I.D.)  
Exterior 1.05" Ht x 3.25" O.D. (26.7 mm Ht x 82.6 mm O.D.)  
Mesh Material: SS 316 / Mesh Size: 35x35 / Opening: 0.017" (0.43mm)



**2" Large Capacity Holder (Material: SS316) - Cat# 8771B**  
Interior 0.7" Ht x 2.0" I.D. (17.8 mm Ht x 50.8 mm I.D.)  
Exterior 1.05" Ht x 2.25" O.D. (26.7 mm Ht x 57.2 mm O.D.)  
Mesh Material: SS 316 / Mesh Size: 35x35 / Opening: 0.017" (0.43mm)



**1" Large Capacity Holder (Material: SS316) - Cat# 8771A**  
Interior 0.7" Ht x 1.0" I.D. (17.8 mm Ht x 25.4 mm I.D.)  
Exterior 1.05" Ht x 1.20" O.D. (26.7 mm Ht x 30.5 mm O.D.)  
Mesh Material: SS 316 / Mesh Size: 35x35 / Opening: 0.017" (0.43mm)

**4" Large Capacity Holder (Cat# 8771D)**



## 3" Large Capacity Holder (Cat# 8771C)




**2" Large Capacity Holder (Cat# 8771B)**



## 1" Large Capacity Holder (Cat# 8771A)



**Opening Large Capacity Holders**



# Large Capacity Holders

---

