

SOTER™ CONDENSER (#8777)



Note: Overall System Set-Up Area Requirement with optional SOTER™ Condenser is 36”(W) x 30”(D) 92cm(W) x 76cm(D)



8777 SOTER™ Condenser

SOTER™ Condenser capable of quietly separating alcohol(s) from CO₂ exhaust discharge of your SAMDRI® chamber. Comes complete with 10' static free 1/4" ID exhaust line, insulated exhaust connect from SAMDRI®, and clear drain line tubing. 120V/1A and 220V/0.5A versions are available.



SOTER™ Condenser Shown Connected to Autosamdri®-815